

UNCLASSIFIED

AD NUMBER

ADB175871

NEW LIMITATION CHANGE

TO

Approved for public release, distribution
unlimited

FROM

Distribution authorized to U.S. Gov't.
agencies only; Test and Evaluation; 28 JUL
1993. Other requests shall be referred to
Army Medical Research and Development
Command, ATTN: SGRD-RMI-S, Fort Detrick,
MD 21702-5012.

AUTHORITY

AMRMC MCMR-RMI-S [70-1y] memo. dtd 28 Aug
1995

THIS PAGE IS UNCLASSIFIED



DEPARTMENT OF THE ARMY
U.S. ARMY MEDICAL RESEARCH AND MATERIEL COMMAND
FORT DETRICK, MARYLAND 21702-5012



REPLY TO
ATTENTION OF:
MCMR-RMI-S (70-1y)

28 Aug 95

MEMORANDUM FOR Administrator, Defense Technical Information Center, ATTN: DTIC-HDS/William Bush, Cameron Station, Building 5, Alexandria, VA 22304-6145

SUBJECT: Request Change in Distribution Statement

1. The U.S. Army Medical Research and Materiel Command has reexamined the need for the limited distribution statement on technical report for Military Interdepartmental Purchase Request Number 92MM2556. Request the limited distribution statement for AD Number ADB175871 be changed to "Approved for ~~public release~~; distribution unlimited." A copy of this report should be released to the National Technical Information Service.
2. Point of contact for this request is Ms. Judy Pawlus, DSN 343-7322.

per Raylie Rinehart
GARY R. GILBERT
COL, MS
Deputy Chief of Staff
for Information Management

REPRODUCTION QUALITY NOTICE

This document is the best quality available. The copy furnished to DTIC contained pages that may have the following quality problems:

- **Pages smaller or larger than normal.**
- **Pages with background color or light colored printing.**
- **Pages with small type or poor printing; and or**
- **Pages with continuous tone material or color photographs.**

Due to various output media available these conditions may or may not cause poor legibility in the microfiche or hardcopy output you receive.



If this block is checked, the copy furnished to DTIC contained pages with color printing, that when reproduced in Black and White, may change detail of the original copy.

AD-B175 871

AD _____

MIPR NO: 92MM2556

TITLE: PHYSIOLOGICAL STUDY OF 4 MODIFIED SELF-CONTAINED VENTILATORS

SUBTITLE: Evaluation of the Impact 780 Self-Contained Ventilator

PRINCIPAL INVESTIGATOR: Peter H. Abbrecht, M.D., Ph.D.

CONTRACTING ORGANIZATION: Uniformed Services University
of the Health Sciences
4301 Jones Bridge Road
Bethesda, Maryland 20814-4799

REPORT DATE: July 28, 1993

TYPE OF REPORT: Final Report

PREPARED FOR: U.S. Army Medical Research and Development Command, Fort Detrick Frederick, Maryland 21702-5012

DTIC
SELECTED
SEP 01 1993
S B D

DISTRIBUTION STATEMENT: Distribution authorized to U.S. Government agencies only, Test and Evaluation, July 28, 1993. Other requests for this document shall be referred to Commander, U.S. Army Medical Research and Development Command, SGRD-RMI-S, Fort Detrick, Frederick, Maryland 21702-5012.

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorized documents.

93-20342

93 8 31 062

REPORT DOCUMENTATION PAGE

GAR 1.1.1

1. AGENCY USE ONLY (Leave Blank)			2. REPORT DATE	3. REPORT TYPE AND DATES COVERED
			28 July 1993	Final Report (4/1/92 - 12/31/92)
4. TITLE AND SUBTITLE PHYSIOLOGICAL STUDY OF 4 MODIFIED SELF-CONTAINED VENTILATOR SUBTITLE: Evaluation of the Impact 780 Self-Contained Ventilator			5. FUNDING NUMBERS MIPR No. 92MM2556	
6. AUTHOR(S) Peter H. Abbrecht, M.D., Ph.D.				
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Uniformed Services University of the Health Sciences 4301 Jones Bridge Road Bethesda, Maryland 20814-4799			8. PERFORMING ORGANIZATION REPORT NUMBER Project No. G176ER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES) U.S. Army Medical Research & Development Command Fort Detrick Frederick, Maryland 21702-5012			10. SPONSORING/MONITORING AGENCY REPORT NUMBER	
11. SUPPLEMENTARY NOTES				
12a. DISTRIBUTION/AVAILABILITY STATEMENT Distribution authorized to U.S. Government agencies only, Test and Evaluation, July 28, 1993. Other Requests for this document shall be referred to Commander, USAMRDC, Fort Detrick, Frederick, Maryland 21702-5012. <i>ATA:SGRD-RMI-S</i>			12b. DISTRIBUTION CODE	
13. ABSTRACT (Maximum 200 words) We evaluated the performance of four Impact model 780 self-contained ventilators before and after vibration testing. Units were tested using the internal compressor and external compressed air powered ventilation modes. We measured impedances to air flow of the input and output ports, the spontaneous inspiration circuit, and the C-2 air filter, and determined the ventilatory parameters including minute ventilation, tidal volume, mean delivered flow rate, insufflation time, exhalation time, and duty cycle. We also determined battery life under the different operating modes. There was relatively little difference in the ventilatory performance of the four different units. The only variable that changed significantly in response to vibration was insufflation time, with a corresponding change in duty cycle. In all four units, the exhalation ports leaked outward during insufflation. All filter connectors leaked unfiltered air inward during spontaneous inspiration and during powered ventilation with the compressor. The average resistance of the four units to inspiratory airflow during spontaneous breathing was 26.5 cm H ₂ O/(L/sec), a value that would probably make spontaneous inspiration through the ventilator impossible for a person with impaired pulmonary function.				
14. SUBJECT TERMS Mechanical Ventilation; Ventilator; Resuscitation Device; R&D V			15. NUMBER OF PAGES	
			16. PRICE CODE	
17. SECURITY CLASSIFICATION OF REPORT Unclassified	18. SECURITY CLASSIFICATION OF THIS PAGE Unclassified	19. SECURITY CLASSIFICATION OF ABSTRACT Unclassified	20. LIMITATION OF ABSTRACT Limited	

GAR 1.1.1

Standard Form 144-1
Prescribed by AFM 10-75
1987-88

Project No: G176ER

**Evaluation of the Impact 780
Self Contained Ventilator**

Peter H. Abbrecht, M.D., Ph.D.

**Uniformed Services University of the
Health Sciences**

**4301 Jones Bridge Road
Bethesda, MD 20814-4799**

APPENDIX - DOCUMENTATION PAGE

is keeping and maintaining the data needed and, consulting and reviewing the content of information, send comments regarding this burden estimate to the Office of Information and Privacy, including suggestions for reducing the burden, to Washington Headquarters Services, Directorate for Information Operations and Resources, 411 Jefferson Davis Highway, Suite 410, Arlington, VA 22204-1202, and to the Office of Management and Budget, Washington, DC 20585.

1. AGENCY USE ONLY (Leave Blank)		2. REPORT DATE	3. REPORT TYPE AND DATES COVERED		
		7-28-1993	Final 4-01-1992 12-31-1993		
4. TITLE AND SUBTITLE				5. FUNDING NUMBERS	
Evaluation of the Impact 780 Self-Contained Ventilator				92MM2556	
6. AUTHOR(S)					
Peter H. Abbrecht, M.D., Ph.D.					
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES)				8. PERFORMING ORGANIZATION REPORT NUMBER	
Uniformed Services University of the Health Sciences 4301 Jones Bridge Rd Bethesda, MD 20814-4799					
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSORING/MONITORING AGENCY REPORT NUMBER	
U.S. Army Medical Materiel Development Activity Ft. Detrick Frederick, MD 21702-5009					
11. SUPPLEMENTARY NOTES					
12a. DISTRIBUTION/AVAILABILITY STATEMENT				12b. DISTRIBUTION CODE	
13. ABSTRACT (Maximum 200 words)					
<p>We evaluated the performance of four Impact model 780 self-contained ventilators before and after vibration testing. Units were tested using the internal compressor and external compressed air powered ventilation modes. We measured impedances to air flow of the input and output ports, the spontaneous inspiration circuit, and the C-2 air filter, and determined the ventilatory parameters including minute ventilation, tidal volume, mean delivered flow rate, insufflation time, exhalation time, and duty cycle. We also determined battery life under the different operating modes. There was relatively little difference in the ventilatory performance of the four different units. The only variable that changed significantly in response to vibration was insufflation time, with a corresponding change in duty cycle. In all four units, the exhalation ports leaked outward during insufflation. All filter connectors leaked unfiltered air inward during spontaneous inspiration and during powered ventilation with the compressor. The average resistance of the four units to inspiratory airflow during spontaneous breathing was 26.5 cm H₂O/(L/sec), a value that would probably make spontaneous inspiration through the ventilator impossible for a person with impaired pulmonary function.</p>					
14. SUBJECT TERMS				15. NUMBER OF PAGES	
Mechanical Ventilation					
				16. PRICE CODE	
17. SECURITY CLASSIFICATION OF REPORT		18. SECURITY CLASSIFICATION OF THIS PAGE		19. SECURITY CLASSIFICATION OF ABSTRACT	
Unclassified		Unclassified		Unclassified	
20. EMBODIMENT CLASSIFICATION				GAR	

The Impact 780 Self-Contained Ventilator consists of an electrically powered compressor, non-adjustable control circuitry, internal battery, external electrical power converter, air filter, hoses, and a multi-valve assembly to join the patient to the ventilator. The internal compressor draws in filtered ambient air and delivers it to the patient. High pressure fittings and a hose are provided to use compressed gas (50 to 80 psi) as an alternate form of powered ventilation. Electrical power is required for either type of powered ventilation, and can be supplied by the external power converter and battery, the external power converter alone, or the battery alone. The patient can always inspire filtered air in addition to, or in place of, either type of powered ventilation. Patient spontaneous inspiration does not require electrical power.

Four ventilators were tested at USUHS to determine the values of the mission critical clinical parameters, both before and after vibration challenge (IAW MIL-STD-810 Vibration, Method 514.3, Procedure X, XI), with three different sources of electrical power. The ventilators were evaluated using the internal compressor and external compressed gas powered ventilation modes. All compressed-gas ventilation testings were performed with compressed air. We also determined the flow/pressure relationships and compliance that would be experienced by a patient breathing spontaneously through the ventilator.

1. General Observations

The exhalation ports on all four patient valves (connector fitting between the patient and the ventilator tubing) leaked outward during powered ventilation insufflation. All filter connectors leaked unfiltered air inward during spontaneous inspiration and during powered ventilation with the internal compressor. In the latter case, from 111 to 765 ml/min of unfiltered air leaked past the filter into the system and was delivered to the patient valve. This leakage around the filter would also expose the patient to unfiltered air during spontaneous inspiration. We were able to eliminate the filter/connector leaks with gaskets. One of the units also had an outward leak inside the ventilator on the positive pressure side of the internal compressor circuit (419 ml/min) which reduced the amount of ventilation delivered to the patient.

The controls and hose fittings on the ventilator were so close together that they were difficult to use. One fuze blew inside a ventilator and another in a power converter. Just to validate the status of the fuses required the use of tools to open the units resulting in exposure of their delicate internal components. The cause(s) of the blown fuses was never discovered. Once the fuses were replaced, there were no further occurrences of this problem.

The external compressed-gas circuit has a low pressure electronic alarm and warning light and a high pressure whistle. The low pressure alarm was active at an average value of 45 psi or less before the vibration challenge, and 47 psi after. The high pressure alarm turned on at an average of 89 psi before the vibration challenge, and 95 psi after. Once activated the high pressure alarm then turned off when the pressure was reduced to an average of 76 psi or less before the vibration challenge, and 84 psi after.

When operating with battery power, battery life was shorter in the internal compressor mode than when using external compressed gas. In the internal compressor mode starting with a fully charged battery, the average time to the first, intermittent low

Dist 3	Avail and/or Special
82	8000

battery alarm buzzer and light was 5.6 hours before the vibration challenge, and 4.8 hours after. The average time to continuous low battery alarm was 6.6 hours before the vibration challenge, and 6.8 hours after. Compressor action stopped completely after an average 7.2 hours operation before the vibration challenge, and 7.9 hours after.

The patient valve has a pressure relief valve to prevent exposing the patient to excessive pressures. The relief valve opened at an average pressure of 31 cm H₂O before the vibration challenge, and 28 cm H₂O after. Escaping flow through the relief valve sounded an alarm at an average pressure of 36 cm H₂O before the vibration challenge, and 35 cm H₂O after. For pressures between the opening level and the alarm level, gas was vented without any audible indication.

During powered ventilation, the maximum pressure that the patient might be exposed to is determined by the amount of flow delivered to the patient valve and by the characteristics of the pressure relief valve. With the internal compressor powered by the external power converter, the unit average peak pressure inside the patient valve was 46 cm H₂O before the vibration challenge and 43 cm H₂O after. With the internal compressor powered by just the internal battery, the average maximum pressure inside the patient valve was 44 cm H₂O before the vibration challenge and 41 cm H₂O after. Using 50 psi external compressed gas, the average maximum patient valve internal pressure was 42 cm H₂O before the vibration challenge and 40 cm H₂O after, and was only slightly influenced by choice of electrical power source.

2. Ventilator Impedance

The amount of ventilation delivered to the patient during powered ventilation depends on the ventilator circuit resistance, the patient's pulmonary resistance, and the amount of power (electrical or gas pressure) driving the ventilator. The following are average resistances at 1 L/sec flow rate measured before and after the vibration challenge.

Vibration Challenge State	<u>Resistance, cm H₂O/(L/sec) ± SD</u>	
	Before	After
Patient Valve Components		
Powered Ventilation Input Port	231 ± 26	264 ± 37
Spontaneous Inspiration Input Port	6.8 ± 0.3	6.3 ± 0.4
Exhalation Output Port, without PEEP	7.5 ± 1.3	9.2 ± 2.5
Exhalation Output Port, with PEEP	19.5 ± 1.9	20.2 ± 2.0
Pressure Relief Valve	72.7 ± 2.6	79.3 ± 5.8
Entire Spontaneous Inspiration Circuit	26.3 ± 1.6	26.7 ± 1.8
Spontaneous Inspiration Hose	6.9 ± 0.5	7.4 ± 0.2
C-2 Air Filter	2.4 ± 0.1	2.4 ± 0.1

The average resistance (26.5 cm H₂O/(L/sec)) of the four units to inspiratory airflow during spontaneous breathing would probably make spontaneous inspiration through the ventilator impossible for a person with impaired pulmonary function.

To test the effects of increased patient impedance on ventilator performance, the systems were tested with seven different resistances at the patient interface port during powered ventilation using both the internal compressor and external compressed air. The following are the average tidal volumes during different powered ventilation modes for two extreme patient resistances.

Added Resistance, cm H ₂ O/(L/sec)	<u>Tidal Volumes, L ± SD</u>	
	0	50
Internal Compressor, Internal Battery and External Power Converter	0.96 ± 0.06	0.82 ± 0.03
Internal Compressor, Internal Battery	0.82 ± 0.04	0.74 ± 0.03
External Compressed Air at 50 psi Internal Battery and External Power Converter	0.77 ± 0.02	0.72 ± 0.02

3. Ventilator Capability

Compressed-Gas Ventilation

The effect of different inlet pressures on tidal volumes during compressed-air ventilation was determined. The following are the average tidal volumes for supplied inlet air pressures spanning the inlet pressure alarm limits.

	<u>Tidal Volumes, L</u>			
Inlet Pressure, psi	50	60	70	80
Tidal Volumes (liters)	0.76	0.90	1.03	1.18

Spontaneous-Inspiration Compliance

The spontaneous inspiration circuit average static compliance was almost constant for a range of vacuums up to -25 cm H₂O. The average compliance was 0.41 ml/cm H₂O before the vibration challenge and 0.46 ml/cm H₂O after.

Powered Ventilation Performance

The following tables summarize the powered ventilation results for each of the two powered ventilation modes (internal compressor and external compressed air at 50 psi). Each ventilation mode was tested with three difference sources of electrical power (external power converter and battery, external power converter alone, and battery alone). For each ventilation mode and power source, 10 measurements were made on each of the four units. Thus, the mean values are the average of 40 measures (ten replications of each of the four units).

Powered Ventilation Mode: Internal Compressor
Electrical Power Source: External Power Converter and Battery

Before Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	15.48	0.90	0.14
Tidal Volume,	L	0.978	0.054	0.008
Mean Delivered Flow Rate,	L/sec	0.571	0.034	0.005
Frequency,	1/min	15.83	0.32	0.05
Duty Cycle,	%	45.18	1.07	0.17
Insufflation Time,	sec	1.713	0.048	0.008
Expiration Time,	sec	2.078	0.064	0.010

After Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	15.07	0.86	0.14
Tidal Volume,	L	0.959	0.054	0.009
Mean Delivered Flow Rate,	L/sec	0.550	0.029	0.005
Frequency,	1/min	15.72	0.19	0.03
Duty Cycle,	%	45.72	0.68	0.11
Insufflation Time,	sec	1.746	0.031	0.005
Expiration Time,	sec	2.074	0.039	0.006

Powered Ventilation Mode; Internal Compressor
 Electrical Power Source: External Power Converter

Before Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	15.15	0.92	0.15
Tidal Volume,	L	0.960	0.056	0.009
Mean Delivered Flow Rate,	L/sec	0.561	0.037	0.006
Frequency,	1/min	15.79	0.34	0.05
Duty Cycle,	%	45.03	1.00	0.16
Insufflation Time,	sec	1.712	0.058	0.009
Expiration Time,	sec	2.090	0.052	0.008

After Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	14.71	0.74	0.12
Tidal Volume,	L	0.939	0.048	0.008
Mean Delivered Flow Rate,	L/sec	0.537	0.026	0.004
Frequency,	1/min	15.67	0.21	0.03
Duty Cycle,	%	45.74	0.72	0.12
Insufflation Time,	sec	1.752	0.033	0.005
Expiration Time,	sec	2.079	0.042	0.007

Powered Ventilation Mode: Internal Compressor
Electrical Power Source: Battery

Before Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.89	0.73	0.12
Tidal Volume,	L	0.824	0.046	0.007
Mean Delivered Flow Rate,	L/sec	0.483	0.024	0.004
Frequency,	1/min	15.64	0.24	0.04
Duty Cycle,	%	44.55	1.13	0.18
Insufflation Time,	sec	1.710	0.045	0.007
Expiration Time,	sec	2.128	0.060	0.009

After Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.78	0.56	0.09
Tidal Volume,	L	0.817	0.036	0.006
Mean Delivered Flow Rate,	L/sec	0.473	0.019	0.003
Frequency,	1/min	15.64	0.23	0.04
Duty Cycle,	%	45.02	0.85	0.13
Insufflation Time,	sec	1.727	0.041	0.006
Expiration Time,	sec	2.110	0.046	0.007

Powered Ventilation Mode: External Compressed-Air (50 psi at inlet during flow)
 Electrical Power Source: External Power Converter and Battery

Before Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.21	0.49	0.08
Tidal Volume,	L	0.767	0.015	0.002
Mean Delivered Flow Rate,	L/sec	0.494	0.016	0.002
Frequency,	1/min	15.91	0.44	0.07
Duty Cycle,	%	41.22	0.79	0.13
Insufflation Time,	sec	1.555	0.034	0.005
Expiration Time,	sec	2.220	0.081	0.013

After Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.18	0.37	0.06
Tidal Volume,	L	0.775	0.021	0.003
Mean Delivered Flow Rate,	L/sec	0.495	0.018	0.003
Frequency,	1/min	15.72	0.20	0.03
Duty Cycle,	%	40.97	0.64	0.10
Insufflation Time,	sec	1.565	0.030	0.005
Expiration Time,	sec	2.255	0.039	0.006

Powered Ventilation Mode: External Compressed-Air (50 psi at inlet during flow)
Electrical Power Source: External Power Converter

Before Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.12	0.39	0.06
Tidal Volume,	L	0.767	0.015	0.002
Mean Delivered Flow Rate,	L/sec	0.493	0.017	0.003
Frequency,	1/min	15.81	0.31	0.05
Duty Cycle,	%	40.95	0.77	0.12
Insufflation Time,	sec	1.555	0.042	0.007
Expiration Time,	sec	2.243	0.052	0.008

After Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.21	0.35	0.06
Tidal Volume,	L	0.778	0.020	0.003
Mean Delivered Flow Rate,	L/sec	0.497	0.017	0.003
Frequency,	1/min	15.68	0.20	0.03
Duty Cycle,	%	40.94	0.58	0.09
Insufflation Time,	sec	1.567	0.031	0.005
Expiration Time,	sec	2.261	0.035	0.006

Powered Ventilation Mode: External Compressed-Air (50 psi at inlet during flow)
Electrical Power Source: Battery

Before Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.04	0.37	0.06
Tidal Volume,	L	0.767	0.0014	0.002
Mean Delivered Flow Rate,	L/sec	0.491	0.014	0.002
Frequency,	1/min	15.69	0.31	0.05
Duty Cycle,	%	40.91	0.93	0.15
Insufflation Time,	sec	1.566	0.039	0.006
Expiration Time,	sec	2.262	0.064	0.010

After Vibration Challenge

<u>Ventilation Variables</u>	<u>(units)</u>	<u>Mean</u>	<u>Std Dev</u>	<u>Std Error</u>
Minute Ventilation,	L/min	12.24	0.38	0.06
Tidal Volume,	L	0.781	0.022	0.003
Mean Delivered Flow Rate,	L/sec	0.498	0.017	0.003
Frequency,	1/min	15.67	0.22	0.04
Duty Cycle,	%	40.97	0.53	0.08
Insufflation Time,	sec	1.568	0.029	0.005
Expiration Time,	sec	2.261	0.038	0.006

The compressor motor used during internal compressor ventilation requires much more power from the battery than does the solenoid during external compressed-gas ventilation. The internal compressor produces a higher flow rate when supplied by the external power converter than when supplied by just the internal battery. Therefore, the internal compressor mode has larger variability for the minute ventilation, tidal volume, and mean delivered flow rate as a function of electrical power source than does the external compressed-gas mode.

The p values of the main effects from the analyses of variance are also listed. The analysis assessed the main effects and their interactions of three independent factors. These factors were: the units, the vibration challenge (VIB), and the sources of electrical power. There were 4 units of the same design. There were two VIB states: before and after. There were three different sources of electrical power: external electrical power converter and battery, external electrical power converter alone, and battery alone. The ten replicated measures were treated as nested terms. We chose p values less than or equal to 0.05 as statistically significant, and they are indicated in bold print.

The following is a listing of p values of the main effects from an analysis of variance for each dependent ventilation variable and each ventilation mode. Each analysis included ten replications for each repeated measure of four units, before and after the vibration challenge, using the three different electric power sources.

	p Values					
	Internal Compressor			External Compressed Air		
	Unit	VIB	Power	Unit	VIB	Power
Minute Ventilation,	0.0001	0.1133	0.0001	0.0001	0.6748	0.0027
Tidal Volume,	0.0001	0.2068	0.0001	0.0001	0.3516	0.1258
Mean Delivered Flow Rate,	0.0001	0.0539	0.0001	0.0001	0.5942	0.5605
Frequency,	0.0001	0.1331	0.0008	0.0001	0.1676	0.0017
Duty Cycle,	0.0001	0.0367	0.0045	0.0001	0.5780	0.1768
Insufflation Time,	0.0001	0.0100	0.0748	0.0001	0.0632	0.1183
Expiration Time,	0.0001	0.3571	0.0007	0.0001	0.2651	0.0045

While there were many statistically significant differences for the various ventilation variables, the differences may not be functionally significant especially because the variation in ventilation requirements of different sized patients may markedly exceed the small differences in ventilator performance observed in these tests.